Twister venturi

No. 1824-0000 / 1824-1000

Bedienungsanleitung
Instruction manual • Mode d'emploi
Istruzioni d'uso • Instrucciones para el servicio
Инструкция по эксплуатации
操作说明书 • 取扱説明書 • 사용 설명서

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ENGLISH

1. Introduction

We are pleased with your decision to purchase a *Twister venturi* vacuum mixer.

Please read the following operating instructions carefully and observe the information they contain in order to ensure a long and trouble-free service life.

2. Application Area

The *Twister venturi* vacuum mixer is designed solely for homogeneous, bubble-free mixing of dental casting and modelling materials such as plasters, investments, and silicones.

2.1 Ambient Conditions

(in accordance with DIN EN 61010)

The unit may only be operated:

- Indoors.
- Up to an altitude of 2,000 m above sea level.
- At an ambient temperature range between 5 - 40°C [41 - 104°F *].
- At a maximum relative humidity of 80% at 31°C [87.8°F], dropping in a linear manner to 50% relative humidity at 40°C [104°F *].
- With mains power where the voltage fluctuations do not exceed 10% of the nominal value.
- Under over-voltage category II conditions.
- Under contamination level 2 conditions.

*) Between 5 – 30°C [41 – 86°F], the unit can be operated at a relative humidity of up to 80%. At temperatures between 31 - 40°C [87.8 - 104°F], the humidity must decrease proportionally in order to ensure operational readiness (e.g., at 35°C [95°F] = 65% humidity; at 40°C [104°F] = 50% humidity). The unit may not be operated at temperatures above 40°C [104°F].

3. Safety Information

3.1 Symbolology

The following symbols are employed in these instructions and on the unit itself:

Danger
This indicates a direct risk of injury. Please observe the operating instructions.

Electrical current
This indicates a hazard due to electrical current.

Attention
Failure to observe the associated information can result in damage to the unit.

Note
This provides the operator with useful information to make working with the unit easier.

Only intended for indoor use.

Before opening the unit, disconnect it from the mains power supply by unplugging the power cord from the wall outlet.

Observe the operating instructions.

3.2 Hazard Information

- The *Twister venturi* is an electric unit and, as such, carries with it a certain inherent potential hazard.
- The unit may only be operated if it is equipped with the proper wall plug and power cord for the mains power system in country of use. Any alterations which may be required to ensure this may only be performed by an electrotechnical specialist.
- The unit may only be operated if the information on the nameplate conforms with the specifications of your local mains power supply.
- The power plug must be easily accessible.
- Prior to use, all electrical and pneumatic lines must be inspected. Units exhibiting damaged lines may not be operated.
- When mixing investment material, please observe the manufacturer's material safety data sheets (hazardous dust) and always wear appropriate protective gear.
- Disconnect the unit from the mains power and compressed air supplies before beginning any maintenance tasks.
- Manipulating the automatic bowl coupler and the intake opening can result in damage to the unit and personal injury.
- The unit may only be used to mix plaster, investment, and silicon.
- Never mix flammable or explosive materials.
- Never connect the blunger without the mixing bowl.
- The use of non-original Renfert blungers and mixing bowls represents a potential injury hazard.

3.3 Liability Exclusion

Renfert shall be absolved from all claims for damages or warranty if:

- The product is employed for any purposes other than those cited in the operating instructions;
- The product is altered in any way other than those alterations described in the operating instructions;
- The product is repaired by other than an authorized facility or if any but Renfert OEM parts are employed;
- The product continues to be employed, despite obvious safety faults or damage;
- The product is subjected to mechanical impacts or is dropped.

4. Installation

4.1 Wall Mounting

(with the included mounting set)

Please have the following on hand:

- Pencil;
- Drilling template;
- Philips head screwdriver;
- Power drill;
- 8 mm [0.3152 inch] drill bit, suitable for penetrating the intended wall material.

! **Inspect the wall on which you intend mounting the unit to ensure that it is capable of bearing the unit's weight!**

1. Determine the optimum working height (Figure 1).
2. Align the drilling template and mark the drill holes (Figure 2).
3. Drill the required holes (using the 8 mm [0.3152 inch] drill bit >>> to a depth of at least 55 mm [2.167 inches]) + insert the dowels (Figure 3).
4. Screw the top screws into the dowels first >>> allowing the screw heads to protrude 9-10 mm [0.3546 - 0.394 inches] (Figure 4).
5. Mount the bottom mounting bracket (Figure 5).
6. Hang the unit on the wall (Figure 6).
7. Secure the unit with the knurled nut >>> Securely tighten the knurled nut (Figure 7).

4.2 Benchtop Units

The *Twister venturi* can be operated as a benchtop unit when a stand is employed.

1. Position the stand on a level surface (Figure 8).
2. Hang the unit in the stand (Figure 9).
3. Secure the unit with the knurled nut >>> Securely tighten the knurled nut (Figure 10).

The stand is not part of the standard delivery and must be ordered as an optional accessory (refer to the „Accessories“ section).

5. Commissioning / Operation

Observe the pressure range information in the technical specifications. A higher inlet pressure will result in a lower vacuum build-up.

For optimal operation as well as to prevent damage, the compressed air should be clean, dry, and oil-free.

Where the compressed air supply produces an inlet pressure that is too high, or where the compressed air is moist, contains oil, or is contaminated, the *Twister venturi* should only be operated with an upstream filter pressure regulator (refer to the „Accessories“ section).

5.1 Operating Elements

(Figure 11 + 12)

- A Vacuum display
- B Vacuum regulator
- C Mixing time display
- D Speed setting
- E ON/OFF switch
- F Control knob (mixing time adjustment, start / stop / venting)
- G Power cord
- H Compressed air connection with inlet filter
- J Vacuum filter
- K Bowl detection pin
- L Silencer
- M Motor fuse
- N Unit fuse

5.2 Commissioning

1. Connect the unit to the compressed air supply (Figure 13a/b). A slight resistance which must be overcome is felt when connecting the compressed air hose to the unit (Figure 13a). Select the proper coupling from the set provided and install it on the other hose end to connect it to the compressed air system (Figure 13b).
2. Connect the unit to the wall outlet (Figure 14).
3. Switch the unit on (Figure 15).
4. Adjust the desired speed >>> 100 - 400 rpm (Figure 16).
5. Adjust the desired mixing time >>> 0 - 5 min. (Figure 17).

Your *Twister venturi* is now operational.

5.3 Mixing Process

When mixing investment material, please always observe the manufacturer's material safety data sheets!

Wear appropriate protective gear as required!

! Do not fill the mixing bowl beyond the maximum fill mark! The maximum fill mark applies to powders and liquids in their un-mixed states.

i Attempting to mix small quantities in large bowls will result in inadequate mixing.

TIP:

Only mix the same type of material in a given bowl. Residue from previous mixing processes can have a negative influence (e.g., silicon fails to harden, etc.).

We recommend a separate bowl for each type of material (plaster, investment, silicon).

1. Select an appropriate size mixing bowl.
2. Using a spatula, lightly mix the materials together (Figure 18).
3. Install the appropriate blunger (Figure 19).

! Make sure the bowl rim (Figure 20a) and cover rim (Figure 20b) are clean!

4. Couple the mixing bowl to the unit (Figure 21). Do not release the bowl until it stays in place by itself.
>>> Indicator outside the red range on the scale (Figure 22a)!

i The motor should start automatically and run briefly when the bowl is inserted. The unit is equipped with an automatic coupling aid that locks the blunger to the motor shaft when the bowl is inserted.

5. Vacuum generation starts automatically.
6. Wait until sufficient vacuum has been created.
7. Start the mixing process
>>> Press the control knob (Figure 22b).
The unit will beep when mixing has been completed.
8. Remove the mixing bowl
>>> Press the control knob (Figure 22b). The bowl will be released in a few seconds.

! Hold the bowl!

5.4 Interrupting the Mixing Process

1. Stop the mixing process
>>> Press the control knob once (Figure 22b).
The unit will beep.
2. Remove the mixing bowl
>>> Press the control knob again (Figure 22b).
The bowl will be released in a few seconds.

! Hold the bowl!

In case of a power failure or if the unit is switched off, the bowl will be vented and be uncoupled.

When power is restored or the unit is again switched on, the mixing motor WILL NOT automatically restart itself.

5.5 Setting the Mixing Time

(Figure 17)

The mixing time set at the control knob (F) is shown on the display (C). The mixing time can be adjusted in 5-second intervals at any time – even during the mixing process itself.

- The mixing time set prior to the start of a mixing process is saved.
- If the mixing time is changed during the course of the mixing process, it is not saved and only applies to the current mixing process.

5.6 Setting the Mixing Speed

(Figure 16)

The mixing speed is adjusted with the adjusting knob (D). The speed can be adjusted at any time – even during the mixing process itself.

5.7 Setting the Vacuum

(Figure 23)

You have the option of individually adjusting the intensity of the vacuum by using the adjusting knob (B).

- + Higher vacuum (greater pressure difference between the pressure in the bowl and the ambient air pressure).
- Lower vacuum (less greater pressure difference between the pressure in the bowl and the ambient air pressure).

Vacuum intensity can also be changed during the mixing process itself.

TIP:

Tests with various investment materials have shown that the best mixing results (smooth, homogeneous casting results) are obtained at the maximum vacuum setting.

As a rule, this also applies to plasters.

In individual instances, independent tests have shown that – under extremely high vacuum – the partial pressure in the mixing bowl can sink to a point where bubbles due to boiling can form in certain plasters (e.g., those containing volatile oils).

Should you observe this phenomenon in the plasters you employ, we recommend first drawing the air out of the mixture for 20 seconds at maximum vacuum pressure, then reducing the vacuum setting to 80-90% during subsequent mixing.

6. Cleaning / Maintenance

Always unplug the unit from the wall outlet before beginning any cleaning or maintenance tasks.

The *Twister venturi* mixing unit is practically maintenance-free.

Nonetheless, the inlet filter and the seals should be cleaned regularly as part of a preventative maintenance programme.

- ! Do not use cleansers containing alcohol or other flammable substances to clean the mixing bowl. Instead, wipe bowls with a damp cloth.

6.1 Cleaning the Inlet Filter

- ! The vacuum built-up is only warranted with cleaned suction filter.

The inlet filter (J) should be inspected and cleaned at least once per month.

The inlet filter should also be cleaned if vacuum fails to build up or is slow in building up.

i The inlet filter is particularly liable to become dirty after mixing silicones or if mixing bowls are completely filled. In this case, you should clean the filter before proceeding with further work.

1. Switch the unit off.
2. Unscrew the inlet filter (Figure 24) using a 17mm spanner.
3. Clean the filter in an ultrasonic bath (we recommend using GO-2011 plaster solvent, Art. No. 2011-0000).
4. Reinstall the inlet filter (Figure 24).

- ! Never operate the unit without the inlet filter!

6.2 Replacing the Intake Filter

In order to prevent dirt particles from being carried into the *Twister venturi* by the compressed air, the compressed air coupling is equipped with an intake filter (H). If vacuum fails to build up or is slow in building up this filter may be blocked and will need to be replaced.

! Always disconnect the unit from the compressed air supply before replacing the intake filter (Figure 13b)

The intake filter is a single unit that is removed and reattached in the same way as a compressed air hose.

- Disconnect the unit from the compressed air supply (Figure 13b).
- Press the ring on the intake filter up and pull the compressed air hose off (Figure 25a).
- Press the ring on the compressed air coupling up and pull the intake filter off (Figure 25b).
- Connect a new intake filter to the compressed air coupling.
- Connect the compressed air hose to the new intake filter.
- Reconnect the unit to the compressed air supply (Figure 13b).
- Gently pull on the intake filter and the hose to make sure all connections are secure.

i Frequent intake filter blockage is an indication of contaminated compressed air. In such cases we recommend the upstream installation of a filter pressure regulator (refer to the „Accessories“ section).

6.3 Seal Surfaces

The following seal surfaces must always be kept clean in order to ensure correct vacuum built up and the secure retention of the mixing bowl during the mixing process:

- Seal between the unit and the blunger (Figure 26a).
- Seal between the blunger and the bowl (Figure 26b).

TIP:

You should coat the rubber seals on the blunger with Vaseline once every month to ensure optimum vacuum force.

6.4 Replacing the Electric Unit and Motor Fuse

1. Unplug the unit from the wall outlet (Figure 14).
2. Unscrew the fuse holder cover (Figure 12 (M), (N), Figure 27a + b).
3. Remove the fuse and replace it (Figure 28). Refer to the „Technical Specifications“ for the proper fuse ratings.

Never use fuses with a higher current rating!

4. Reinstall and close the fuse holder (Figure 27).

7. Spare Parts

For the order numbers of consumables and spare parts please refer to the spare parts list at the end of this manual.

8. Warranty

Provided the unit is properly used, Renfert **warrants** the *Twister venturi* for a period of **3 years**.

Presentation of the original bill of sale is required for all warranty claims.

Components subject to natural wear (e.g., blungers, mixing bowls, and intake filters) are excluded from this warranty. The warranty is voided in case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in case of independent repairs or repairs by unauthorized personnel; if spare parts from other manufacturers are employed, or; in case of unusual influences or influences not in compliance with the utilization instructions. Warranty service shall not extend the original warranty.

9. Technical Specifications

| | |
|------------------------------|--|
| Mains voltage: | 230 V, 50-60 Hz 120 V, 50-60 Hz |
| Power consumption: | 150 W |
| Device fuse: | T 1.6 A (230 V) T 3.15 A (120 V) |
| Motor fuse: | T 8 A |
| Speed: | 150 to 400 rpm |
| Vacuum, max.: | approx. -880 mbar; corresponds to a bowl pressure of approx. 90 mbar |
| Capacity, max.: | 24 l/min |
| Air consumption, approx.: | 46 l/min at 4.5 bar dynamic inlet pressure |
| Operating pressure, dynamic: | 4.5 to 7 bar |
| Dimensions (WxHxD): | 140 x 325 x 225 mm [5.51 x 12.81 x 8.87 inches] |
| Weight: | 6.5 kg, w./o. bowl |

10. Standard Delivery

- 1 *Twister venturi* vacuum mixing unit
- 1 500 ml bowl, incl. blunger
- 1 Compressed air hose, 2 m
- 1 Pneumatic connection set
- 1 Set of operating instructions
- 1 Spare parts list
- 1 Drilling template
- 1 Fastening set

11. Delivery Versions

- No. 1824-0000 *Twister venturi* 230 V / 50-60 Hz, incl. 500 ml-bowl with blunger
- No. 1824-1000 *Twister venturi* 120 V / 50-60 Hz, incl. 500 ml-bowl with blunger

12. Accessories

- No. 2929-0000 Filter pressure regulator
- No. 1821-0101 Stand for benchtop units;
WxHxD: 230 x 680 x 290 mm
[9.062 x 26.772 x 11.426 inches]
- No. 1821-0200 Mixing spatula
- No. 1820-6500 Bowl incl. blunger, 65 ml
- No. 1820-6510 Blunger, 65 ml
- No. 1820-6520 Bowl, 65 ml
- No. 1820-0200 Bowl, incl. blunger, 200 ml
- No. 1820-0210 Blunger, 200 ml
- No. 1820-0220 Bowl, 200 ml
- No. 1820-0500 Bowl, incl. blunger, 500 ml
- No. 1820-0510 Blunger, 500 ml
- No. 1820-0520 Bowl, 500 ml
- No. 1820-0700 Bowl, incl. blunger, 700 ml
- No. 1820-0710 Blunger, 700 ml
- No. 1820-0720 Bowl, 700 ml
- No. 1820-1001 Bowl, incl. blunger, 1000 ml
- No. 1820-1010 Blunger, 1000 ml
- No. 1820-1020 Bowl, 1000 ml

13. Disposal instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste:

For more information regarding proper disposal please apply at your local authorities.

14. Troubleshooting

| Problem | Possible cause | Corrective action |
|--|---|--|
| ON/OFF switch fails to work. | <ul style="list-style-type: none"> • Power not connected. • Device fuse blown. • Damaged power cord. | <ul style="list-style-type: none"> • Check the power supply. • Replace the fuse (refer to Sec. 6.4). • Have the device repaired. |
| Motor fails to start. | <ul style="list-style-type: none"> • Motor fuse blown. • Motor fault. | <ul style="list-style-type: none"> • Replace the fuse (refer to Sec. 6.4). • Have the device repaired. |
| No, or insufficient vacuum (not in the green range on the scale). | <ul style="list-style-type: none"> • Intake filter plugged. • Seal surfaces dirty. • Compressed air hose not connected, leaky or bent. • Venturi nozzle plugged. • Intake filter plugged • Compressor switched off. • Solenoid valve fault. • Silencer plugged. • Vacuum adjusting knob turned out too much. • Compressed air system with insufficient performance. • Working pressure too low. • Air hose too long. • Square are of air hose too small. | <ul style="list-style-type: none"> • Clean the filter (refer to Sec. 6.1). • Clean the seal surfaces (refer to Sec. 6.3). • Check the pipe. • Have the device repaired. • Replace intake filter (refer to Sec. 6.2). • Switch the compressor on. • Have the device repaired. • Replace the silencer. • Adjust vacuum (see pt. 5.7). • Control the compressor's performance (see technical specifications, eventually switch other users off). • Increase dynamic working pressure (4,5 to 7 bar / 65,27 to 101,53 psi.). • Maximum length 2 m. • Minimum inner diameter 4 mm. |
| Vacuum release / venting, too slow. | <ul style="list-style-type: none"> • Intake filter plugged. • Silencer plugged. • Solenoid valve fault | <ul style="list-style-type: none"> • Clean the filter (refer to Sec. 6.1). • Clean resp. replace silencer. • Connect a filter unit. • Have the device repaired. |
| Permanent venting noise. | <ul style="list-style-type: none"> • Bowl detection pin sticking. | <ul style="list-style-type: none"> • Have the device repaired. |
| Paddle does not rotate. | <ul style="list-style-type: none"> • Paddle is overloaded, motor fuse (8A) is blown. | <ul style="list-style-type: none"> • Clean the paddle and the bowl. • Replace the fuse (refer to Sec. 6.4). |
| Fluttering noise from the silencer. | <ul style="list-style-type: none"> • Working pressure too low | <ul style="list-style-type: none"> • Increase dynamic working pressure (4,5 to 7 bar / 65,27 to 101,53 psi). |

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Renfert GmbH • Industriegebiet • 78247 Hilzingen/Germany